

LABORATORY REPORT

3–5

Student: _____

Desk No.: _____ Section: _____

Ex. 3 Phase-Contrast Microscopy

A. Questions

1. Which rays (*direct* or *diffracted*) are altered by the phase ring on the phase plate? _____
2. How much phase shift occurs in the light rays that emerge from a transparent object? _____
3. Differentiate:

Bright phase microscope: _____

Dark phase microscope: _____

4. List two items that can be used for observing the concentricity of the annulus and phase ring:
a. _____ b. _____

Ex. 4 Fluorescence Microscopy

A. Questions

1. Differentiate:

Phosphorescence: _____

Fluorescence: _____

2. List three fluorochromes that are used in staining bacteria.

a. _____ b. _____ c. _____

3. What are the two most serious hazards when using mercury vapor arc lamps?

a. _____ b. _____

4. What relationship exists between the wavelength of light and its energy?

Ex. 5 Microscopic Measurements

A. Questions

1. What is the distance between each of the graduations on the stage micrometer? _____ mm
2. Why must the entire calibration procedure be performed for each objective? _____

(See reverse page for more questions pertaining to Exercises 3 and 4.)

Phase-Contrast and Fluorescence Microscopy

B. Multiple Choice Questions for Exercises 3 and 4

Select the answer that best completes the following statements.

1. If direct rays passing through an object are advanced $\frac{1}{4}$ wavelength by the phase ring, the diffracted rays are
 1. in phase with the direct rays.
 2. $\frac{1}{2}$ wavelength out of phase with the direct rays.
 3. $\frac{3}{4}$ wavelength out of phase with direct rays.
 4. in reverse phase with the direct rays.
 5. Both 2 and 4 are correct.
2. The best barrier filter to use with the Schott BG12 exciter filter is the
 1. blue A0702.
 2. Schott OG1.
 3. Wratten G-2.
 4. None of these are correct.
3. The examination of ordinary stained slides is not enhanced with a
 1. brightfield microscope.
 2. phase-contrast microscope.
 3. fluorescence microscope.
 4. Both 2 and 3 are correct.
4. Special immersion oil is required for
 1. brightfield microscopy.
 2. phase-contrast microscopy.
 3. fluorescence microscopy.
 4. Both 2 and 3 are correct.
5. Amplitude summation occurs in phase-contrast optics when both direct and diffracted rays are
 1. in phase.
 2. in reverse phase.
 3. off $\frac{1}{4}$ wavelength.
 4. None of these are correct.
6. The phase-contrast microscope is best suited for observing
 1. living organisms in an uncovered drop on a slide.
 2. stained slides with cover glasses.
 3. living organisms in hanging drop slide preparations.
 4. living organisms on a slide with a cover glass.
7. The wavelength of fluorescent light rays is
 1. always longer than the exciting wavelength.
 2. always shorter than the exciting wavelength.
 3. about the same length as the exciting wavelength.
 4. sometimes shorter and at other times longer in wavelength.
8. The barrier filter in the fluorescence microscope is kept in position to
 1. block all light rays from getting through.
 2. allow fluorescence light rays to pass through.
 3. screen out exciting light rays.
 4. Both 2 and 3 are correct.
9. A phase centering telescope is used to
 1. improve the resolution of the ocular.
 2. increase magnification with the oil immersion objective.
 3. observe the relationship of the annular diaphragm to the phase ring.
 4. None of the above are correct.
10. The visible spectrum of light is between
 1. 200 and 800 nanometers.
 2. 400 and 780 nanometers.
 3. 300 and 800 nanometers.
 4. None of these are correct.

Answers

Multiple Choice

1. _____
2. _____
3. _____
4. _____
5. _____
6. _____
7. _____
8. _____
9. _____
10. _____

LABORATORY REPORT

6

Student: _____

Desk No.: _____ Section: _____

Protozoa, Algae, and Cyanobacteria

A. Tabulation of Observations

In this study of freshwater microorganisms, record your observations in the following tables. The number of organisms to be identified will depend on the availability of time and materials. Your instructor will indicate the number of each type that should be recorded.

Record the genus of each identifiable type. Also, indicate the phylum or division to which the organism belongs. Microorganisms that you are unable to identify should be sketched in the space provided. It is not necessary to draw those that are identified.

PROTOZOA

GENUS	PHYLUM	BOTTLE NO.	SKETCHES OF UNIDENTIFIED

ALGAE

GENUS	DIVISION	BOTTLE NO.	SKETCHES OF UNIDENTIFIED

Protozoa, Algae, and Cyanobacteria

CYANOBACTERIA

GENUS	BOTTLE NO.	SKETCHES OF UNIDENTIFIED

B. General Questions

Record the answers to the following questions in the answer column. It may be necessary to consult your text or library references for one or two of the answers.

- Give the kingdom in which each of the following groups of organisms is found:
 - protozoans
 - algae
 - cyanobacteria
 - bacteria
 - fungi
 - microscopic invertebrates
- Four kingdoms are represented by the organisms in the above question. Name the fifth kingdom.
- What is the most significant characteristic seen in eukaryotes that is lacking in prokaryotes?
- What characteristic in the microscopic invertebrates distinguishes them from protozoans?
- Which protozoan phylum was not found in pond samples because phylum members are all parasitic?
- Indicate whether the following are *present* or *absent* in the algae:
 - cilia
 - flagella
 - chloroplasts
- Indicate whether the following are *present* or *absent* in the protozoans:
 - cilia
 - chloroplasts
 - mitochondria
 - mitosis
- Which photosynthetic pigment is common to all algae and cyanobacteria?
- Name two photosynthetic pigments that are found in the cyanobacteria but not in the algae.
- What photosynthetic pigment is found in bacteria but is lacking in all other photosynthetic organisms?
- What type of movement is exhibited by the diatoms?

Answers

- _____
 - _____
 - _____
 - _____
 - _____
 - _____
- _____
- _____
- _____
- _____
- _____
 - _____
 - _____
- _____
 - _____
 - _____
 - _____
- _____
- _____
 - _____
- _____
- _____
- _____

Protozoa, Algae, and Cyanobacteria

C. Protozoan Characterization

Select the protozoan groups in the right-hand column that have the following characteristics:

- | | |
|------------------------------|------------------|
| 1. move with flagella | 1. Sarcodina |
| 2. move with cilia | 2. Mastigophora |
| 3. move with pseudopodia | 3. Ciliophora |
| 4. have nuclear membranes | 4. Sporozoa |
| 5. lack nuclear membranes | 5. all of above |
| 6. all species are parasitic | 6. none of above |
| 7. produce resistant cysts | |

D. Characterization of Algae and Cyanobacteria

Select the groups in the right-hand column that have the following characteristics:

Pigments

- | | |
|--------------------|----------------------|
| 1. chlorophyll a | 1. Euglenophycophyta |
| 2. chlorophyll b | 2. Chlorophycophyta |
| 3. chlorophyll c | 3. Chrysophycophyta |
| 4. fucoxanthin | 4. Phaeophycophyta |
| 5. c-phycoerythrin | 5. Pyrrophytocypha |
| 6. c-phycoerythrin | 6. Cyanobacteria |
| | 7. all of above |
| | 8. none of above |

Food Storage

7. fats
8. oils
9. starches
10. laminarin
11. leucosin
12. paramylum
13. mannitol

Other Structures

14. pellicle, no cell wall
15. cell walls, box and lid
16. chloroplasts
17. phycobilisomes
18. thylakoids

Answers

Protozoa

1. _____
2. _____
3. _____
4. _____
5. _____
6. _____
7. _____

Algae

1. _____
2. _____
3. _____
4. _____
5. _____
6. _____
7. _____
8. _____
9. _____
10. _____
11. _____
12. _____
13. _____
14. _____
15. _____
16. _____
17. _____
18. _____

Benson: Microbiological Applications Lab Manual, Eighth Edition	Back Matter	Laboratory Reports		© The McGraw-Hill Companies, 2001
--	--------------------	---------------------------	--	--------------------------------------

LABORATORY REPORT

8, 9

Student: _____

Desk No.: _____ Section: _____

Ex. 8 Aseptic Technique

A. Results

1. Were all your transfers successful? _____
2. What evidence do you have that they were successful? _____

3. What evidence do you have that a transfer is unsuccessful? _____

B. Questions:

1. What kinds of organisms are destroyed when your desktop is scrubbed down with a disinfectant? _____

2. Are bacterial endospores destroyed? _____
3. How hot should inoculating loops and needles be heated? _____
4. Why is it necessary to flame the mouth of the tube before and after performing an inoculation? _____

Ex. 9 The Bacteria

A. Questions:

After you have tabulated your results on the back of this sheet, answer the following questions:

1. Using the number of colonies as an indicator, which habitat sampled by the class appears to be the most contaminated one? _____
2. Why do you suppose this habitat contains such a high microbial count? _____

3.
 - a. Were any plates completely lacking in colonies? _____
 - b. Do you think that the habitat sampled was really sterile? _____
 - c. If your answer to *b* is *no*, then how can you account for the lack of growth on the plate? _____

 - d. If your answer to *b* is *yes*, defend it: _____

4. In a few words describe some differences in the macroscopic appearance of bacteria and mold colonies: _____

Bacteria

B. Tabulation

After examining your TSA and blood agar plates, record your results in the following table and on a similar table that your instructor has drawn on the chalkboard. With respect to the plates, we are concerned with a quantitative evaluation of the degree of contamination and differentiation as to whether the organisms are bacteria or molds. Quantify your recording as follows:

0 no growth	+	+	+	51 to 100 colonies
+ 1 to 10 colonies	+	+	+	+
+	+	+	+	over 100 colonies

After shaking the tube of broth to disperse the organisms, look for cloudiness (turbidity). If the broth is clear, no bacterial growth occurred. Record no growth as 0. If tube is turbid, record + in last column.

[illegible]

LABORATORY REPORT

10


Student: _____

Desk No.: _____ Section: _____


The Fungi: Yeasts and Molds

A. Yeast Study

Draw a few representative cells of *Saccharomyces cerevisiae* in the appropriate circles below. Blastospores (buds) and ascospores, if seen, should be shown and labeled.



Prepared Slide



Living Cells

B. Mold Study

In the following table, list the genera of molds identified in this exercise. Under colony description, give the approximate diameter of the colony, its topside color and backside (bottom) color. For microscopic appearance, make a sketch of the organism as it appears on slide preparation.

GENUS	COLONY DESCRIPTION	MICROSCOPIC APPEARANCE (DRAWING)

Fungi: Yeasts and Molds

C. Questions

Record the answers for the following questions in the answer column.

1. The science that is concerned with the study of fungi is called _____.
2. The kingdom to which the fungi belong is _____.
3. Microscopic filaments of molds are called _____.
4. A filamentlike structure formed by a yeast from a chain of blastospores is called a _____.
5. A mass of mold filaments, as observed by the naked eye, is called a _____.
6. Most molds have _____ hyphae (*septate* or *non-septate*).
7. List three kinds of sexual spores that are the basis for classifying the molds.
8. What is the name of the rootlike structure that is seen in *Rhizopus*?
9. What type of hypha is seen in *Mucor* and *Rhizopus*?
10. What kind of asexual spores are seen in *Mucor* and *Rhizopus*?
11. What kind of asexual spores are seen in *Penicillium*?
12. What kind of asexual spores are seen in *Alternaria*?
13. Which subdivision of the Amastigomycota contains individuals that lack sexual spores?
14. What division of Myceteae consists of slime molds?
15. Fungi that exist both as yeasts and molds are said to be _____.

Answers

1. _____
2. _____
3. _____
4. _____
5. _____
6. _____
- 7.a. _____
- b. _____
- c. _____
8. _____
9. _____
10. _____
11. _____
12. _____
13. _____
14. _____
15. _____

LABORATORY REPORT

11–14

Student: _____


Desk No.: _____ Section: _____

Negative Staining, Smear Preparation, Simple and Capsular Staining

A. Negative Staining (Exercise 11)

1. **Drawing:** Make a drawing in the circle at the right of some of the organisms as seen under oil immersion.
2. In addition to nigrosine, what other agent is often used for making negative-stained slides?

3. Other than bacteria, what kinds of microorganisms might one encounter in the mouth?



Oral organisms
(nigrosine)


B. Smear Preparation (Exercise 12)

1. Give two reasons for heating the slide after the smear is air-dried:
 - a. _____
 - b. _____
2. Why is an inoculating needle preferred to a wire loop when making smears from solid media?

C. Simple Staining (Exercise 13)

1. **Drawing:** Draw a few cells of *C. diphtheriae* from the portion of the slide that exhibits metachromatic granules and palisade arrangement.
2. Why are basic dyes more successful on bacteria than acidic dyes?

3. List three basic dyes that are used to stain bacteria: a. _____ ,
b. _____ , c. _____




Corynebacterium diphtheriae

Negative Staining, Smear Preparation, Simple and Capsular Staining

D. Capsular Staining (Exercise 14)

1. List some of the chemical substances that have been identified in bacterial capsules.

2. What relationship is there between capsules and bacterial virulence?



Klebsiella pneumoniae
(capsular stain)

LABORATORY REPORT

15–18


Student: _____

Desk No.: _____ Section: _____


Differential Staining: Gram, Spore, and Acid-Fast

A. Drawings


With colored pencils, draw the various organisms as seen under oil immersion. Extra circles are provided for additional assignments, if needed.




P. aeruginosa & *S. aureus*
(Gram stain)




B. megaterium & *M. B. catarrhalis*
(Gram stain)




M. smegmatis
(Gram stain)




B. megaterium
(Schaeffer-Fulton method)



B. megaterium
(Dorner method)



M. smegmatis & *S. aureus*
(Ziehl-Neelsen method)



M. phlei & *S. aureus*
(Truant method)

OPTIONAL STAINING

Differential Staining: Gram, Spore, and Acid-Fast

B. Completion Questions

1. What color would you expect *S. aureus* to be if the iodine step were omitted in the Gram staining procedure? _____
Explanation: _____

2. What part of the bacterial cell (*cell wall* or *protoplast*) appears to play the most important role in determining whether an organism is gram-positive?

3. Why would methylene blue not work just as well as safranin for counterstaining in the Gram staining procedure? _____
4. Why are endospores so difficult to stain? _____
5. How do the following two genera of spore-formers differ physiologically?
Bacillus: _____
Clostridium: _____
6. How do you differentiate *S. aureus* and *M. B. catarrhalis* from each other on the basis of morphological characteristics? _____

7. Are the acid-fast mycobacteria gram-positive or gram-negative? _____
8. For what two diseases is acid-fast staining of paramount importance?
a. _____ b. _____
9. What advantage does the Ziehl-Neelsen technique have over the Truant method?

10. What advantage does the Truant method have over the Ziehl-Neelsen technique?

11. Why is it desirable to combine *S. aureus* with acid-fast organisms such as *M. smegmatis* when applying an acid-fast staining technique?

LABORATORY REPORT

19

Student: _____

Desk No.: _____ Section: _____

Motility Determination

A. Test Results

1. Which of the two organisms exhibited true motility on the slides?

2. Did the semisolid medium inoculations confirm the results obtained from the slides?

3. Sketch in the appearance of the two tube inoculations:



Micrococcus luteus



Proteus vulgaris

B. Questions

1. How does Brownian movement differ from true motility? _____

2. How do you differentiate water current movement from true motility? _____

3. Make sketches that illustrate each of the following flagellar arrangements:

Monotrichic

Lophotrichic

Amphitrichic

Peritrichic

Benson: Microbiological Applications Lab Manual, Eighth Edition	Back Matter	Laboratory Reports		© The McGraw-Hill Companies, 2001
--	--------------------	---------------------------	--	--------------------------------------

LABORATORY REPORT

20

Student: _____

Desk No.: _____ Section: _____

Culture Media Preparation

1. How do the following types of organisms differ in their carbon needs?

Photoautotrophs: _____

Photoheterotrophs: _____

2. Where do the above two types of organisms get their energy? _____

3. Where do **chemoheterotrophs** get their energy? _____

4. What is a growth factor? _____

5. Give two reasons why agar is such a good ingredient for converting liquid media to solid media.

a. _____

b. _____

6. Differentiate between the following two types of media:

Synthetic medium: _____

Nonsynthetic medium: _____

7. Differentiate between the following two types of media:

Selective medium: _____

Differential medium: _____

8. Briefly, list the steps that you would go through to make up a batch of nutrient agar slants:

Benson: Microbiological Applications Lab Manual, Eighth Edition	Back Matter	Laboratory Reports		© The McGraw-Hill Companies, 2001
--	--------------------	---------------------------	--	--------------------------------------

LABORATORY REPORT

21


Student: _____

Desk No.: _____ Section: _____

Pure Culture Techniques

A. Evaluation of Streak Plate

Show within the circle the distribution of the colonies on your streak plate. To identify the colonies, use red for *Serratia marcescens*, yellow for *Micrococcus luteus*, and purple for *Chromobacterium violaceum*. If time permits, your instructor may inspect your plate and enter a grade where indicated.



Grade _____

B. Evaluation of Pour Plates

Show the distribution of colonies on plates II and III, using only the quadrant section for plate II. If plate III has too many colonies, follow the same procedure. Use colors.




plate II




plate III


Pure Culture Techniques

C. Subculture Evaluation


With colored pencils, sketch the appearance of the growth on the slant diagrams below. Also, draw a few cells of each organism as revealed by Gram staining in the adjacent circle.




Serratia marcescens



Micrococcus luteus
or
Chromobacterium violaceum



Escherichia coli



D. Questions

1. Which method of separating organisms seems to achieve the best separation?

2. Which method requires the greatest skill? _____

3. Do you think you have pure cultures of each organism on the slants? _____

Can you be absolutely sure by studying its microscopic appearance? _____

Explain: _____

4. Give two reasons why the nutrient agar must be cooled to 50° C before inoculating and pouring.

5. Why should a Petri plate be discarded if media is splashed up the side to the top?

6. Give two reasons why it is important to invert plates during incubation:

7. Why is it important not to dig into the agar with the loop? _____

8. Why must the loop be flamed before entering a culture? _____

Why must it be flamed after making an inoculation? _____

LABORATORY REPORT

22


Student: _____

Desk No.: _____ Section: _____


Cultivation of Anaerobes

A. Tube Inoculations

After carefully comparing the appearance of the six cultures belonging to you and your laboratory partner, select the best tube for each organism and sketch its appearance in the tubes below. Indicate under each name the type of medium (FTM or TGYA).



E. coli
()




S. faecalis
()



S. aureus
()



B. subtilis
()



C. sporogenes
()



C. rubrum
()

B. Plate Inoculations

After comparing the growths on the two plates of Brewer's anaerobic agar with the growths in the six tubes, classify each organism as to its oxygen requirements:

Escherichia coli: _____

Bacillus subtilis: _____

Streptococcus faecalis: _____

Clostridium sporogenes: _____

Staphylococcus aureus: _____

Clostridium rubrum: _____

C. Questions

1. What is the function of oxygen at the cellular level? _____

2. Why are facultative organisms able to grow with or without oxygen while aerobes grow only in its presence? _____

3. How do "indifferents" differ from "facultatives"? _____


4. What is the function of the following agents in the media used in this experiment?
 Sodium thioglycollate: _____
 Resazurin: _____
 Agar in FTM: _____
 Agar in TGYA shake: _____

Cultivation of Anaerobes


5. How is oxygen removed from the air in a GasPak anaerobic jar? _____

D. Spore Study


If a spore-stained slide is made of the three spore-formers, draw a few cells of each organism in the spaces provided below:



B. subtilis



C. rubrum



C. sporogenes

LABORATORY REPORT

23

Student: _____

Desk No.: _____ Section: _____

Bacterial Population Counts

A. Quantitative Plating Method

1. Record your plate counts in this table:

DILUTION PLATED	ML PLATED	NUMBER OF COLONIES
1:10,000	1.0	
1:100,000	0.1	
1:1,000,000	1.0	
1:10,000,000	0.1	

2. How many cells per ml were there in the undiluted culture? _____
3. How would you inoculate a plate to get 1:100 dilution? _____

4. How would you inoculate a plate to get 1:10 dilution? _____

5. Give two reasons why it is necessary to shake the water blanks as recommended.
- a. _____
- b. _____


B. Turbidimetric Determinations

1. Record the percent transmittance and optical density values for your dilutions in the following table.

DILUTION	PERCENT TRANSMITTANCE	OPTICAL DENSITY
1:1		
1:2		
1:4		
1:8		
1:16		

Bacterial Population Counts

- Plot the optical densities versus the concentration of organisms. Complete the graph by drawing a line between plot points.



C. Questions

- What is the maximum O.D. that is within the linear portion of the curve? _____
- What is the corrected or true O.D. of the undiluted culture? (Hint: If the O.D. for the 1:2 dilution but not the 1:1 dilution is within the linear portion of the curve, then the O.D. of the 1:1 dilution should not be considered correct. The correct or true O.D. of the undiluted culture in this example could be estimated by multiplying the O.D. of the 1:2 dilution by 2.) _____
- What is the correlation between corrected O.D. and cell number for your culture? _____

- Why is it necessary to perform a plate count in conjunction with the turbidimetry procedure? _____

- If your medium were pale blue instead of amber-colored, as is the case of nutrient broth, would you set the wavelength control knob higher or lower than 686 nanometers? _____

LABORATORY REPORT

24, 25


Student: _____

Desk No.: _____ Section: _____


Ex. 24 Slide Culture: Autotrophs

A. Microscopic Examination

While examining the two slides, move them around to different areas to note the various types of organisms that are present. Draw representative types.



GRAM'S STAIN



LIVING

B. Questions

1. With respect to Gram's stain, which type (*gram-positive* or *gram-negative*) seems to predominate?

2. List as many different kinds of autotrophic protists as you can that can be cultured on this type of slide.

3. Some organisms that grow on this type of slide are chemosynthetic heterotrophs. What would be the source of their food?

Slime Mold Culture

Ex. 25 Slime Mold Culture

A. Observations

1. What happened when the flow of protoplasm on a plasmodium was interrupted by severance with a scalpel?

2. Describe your observations of the crushed spores on the hanging drop slide.

B. Questions

1. List two functions served by fructification (sporangia formation) in *Physarum*.

a. _____

b. _____

2. What is the principal function of the plasmodial stage of *Physarum*?

3. List two characteristics that the Myxobacterales and Gymnomycota have in common.

a. _____

b. _____

4. Postulate as to the evolutionary relationship between the myxobacteria and slime molds.

LABORATORY REPORT

27

Student: _____

Desk No.: _____ Section: _____

Anaerobic Phototrophic Bacteria: Isolation and Culture

A. Observations:

1. **Column Appearance:** Describe in a few words the appearance of the Winogradsky column during the:
First Period:

One Week Later:

Two Weeks Later:

Subsequent Weeks:

2. **Subculture Appearance:** Describe the appearance of your subculture tubes at the end of two weeks:

3. **Microscopic Appearance:** Describe the characteristics of the cells you observed on wet mount slides and gram-stained slides: _____

B. Questions:

1. What roles do the following organisms perform in the Winogradsky column?

Clostridium: _____

Desulfovibrio: _____

2. How do *Chromatium* and *Chlorobium* differ as to where they put the sulfur that they produce? _____

3. Give the equations for photosynthesis in the following:

Algae: _____

***Chromatium* and *Chlorobium*:** _____

Anaerobic Phototrophic Bacteria: Isolation and Culture

C. Tabulation of Results

If different mud samples were used by class members, record your results on the chart below and on a similar chart constructed by your instructor on the chalkboard. Record the presence (+) or absence (–) of the various species identified. In the “Other” column list any other species that you might have encountered.

[illegible]

D. Conclusions from above table:

LABORATORY REPORT

28, 29

Student: _____

Desk No.: _____ Section: _____

Isolation of Phage from Sewage and Flies

A. Plaque Size Increase (for both Exercises 28 and 29)

With a china marking pencil, circle and label three plaques on one of the plates and record their sizes in millimeters at 1-hour intervals.

TIME	PLAQUE SIZE (millimeters)		
	Plaque No. 1	Plaque No. 2	Plaque No. 3
When first seen			
1 hour later			
2 hours later			
3 hours later			

B. Questions (for Exercise 28 only)

1. Were any plaques seen on the negative control plate? _____
2. Do plates 1, 2, and 3 show a progressive increase in number of plaques with increased amount of sewage filtrate? _____
3. Did the phage completely “wipe out” all bacterial growth on any of the plates? _____
If so, which plates? _____

C. Observations (for Exercise 29 only)

Count all the plaques on each plate and record the counts in the following table. If the plaques are very numerous, use a Quebec colony counter and hand counting device. If this exercise was performed as a class project with individual students doing only one or two plates from a common fly-broth filtrate, record all counts on the chalkboard on a table similar to the one below.

Plate Number	1	2	3	4	5	6	7	8	9	10
<i>E. coli</i> (ml)	0.9	0.8	0.7	0.6	0.5	0.4	0.3	0.2	0.1	1.0
Filtrate (ml)	0.1	0.2	0.3	0.4	0.5	0.6	0.7	0.8	0.9	0
Number of plaques										

Isolation of Phage from Sewage and Flies

D. Questions (for Exercise 29 only)

1. Which plate was used as the negative control? _____
2. Were there any plaques on the negative control plate? _____
3. What would be the explanation for the presence of plaques on the negative control plate?

4. Were any plates completely “wiped out” by phage action? _____
If so, which ones? _____

E. Terminology (for both Exercises 28 and 29)

1. Differentiate between the following:
Lysis: _____
Lysogeny: _____
2. Differentiate between the following:
Virulent phage: _____
Temperate phage: _____

LABORATORY REPORT

30

Student: _____

Desk No.: _____ Section: _____

Burst Size Determination: A One-Step Growth Curve

A. Plaque Counts

Record the counts of plaques on each of the plates in the following table. Record the peak number of plaques as the *burst size*. The drop in plaque numbers after a peak results from adsorption of mature phage virions on other bacterial cells and cell debris.

15	25	30	35	40	45	50

Burst size: _____

B. Dilution Interpretation

Answer the following questions to clarify your understanding of the dilutions that occur in this experiment.

- How many cells were present in each milliliter of the original bacterial culture? _____
- How many bacterial cells (total) were dispensed into the ADS tube? _____
- If the bacterial dilution per plate is 1:10,000,000, how many bacterial cells were distributed to each plate?

- How many phage virions were present in 1 ml of the original phage suspension?

- How many phage virions were present in the 0.1 ml of phage suspension that was added to the ADS tube?

- What was the numerical ratio of phage virions to bacterial cells in the ADS tube? _____
What is this ratio called? _____
- How many bacterial cells were placed in the ADS-2 tube? _____
- What effect does dilution have on adsorption? _____

Benson: Microbiological Applications Lab Manual, Eighth Edition	Back Matter	Laboratory Reports		© The McGraw-Hill Companies, 2001
--	--------------------	---------------------------	--	--------------------------------------

LABORATORY REPORT

31–33

Student: _____

Desk No.: _____ Section: _____

Microbial Interrelationships
Ex. 31 Bacterial Commensalism**A. Results**

Indicate the degree of turbidity (none, +, ++, +++) in the following table. With colored pencils, draw the appearance of the gram-stained slides where indicated.

ORGANISMS	TURBIDITY	GRAM STAIN
<i>Staphylococcus aureus</i>		
<i>Clostridium sporogenes</i>		
<i>S. aureus</i> and <i>C. sporogenes</i>		

B. Questions

- Does *C. sporogenes* grow well in nutrient broth? _____
- How does *S. aureus* assist *C. sporogenes* in growth? _____

Ex. 32 Bacterial Synergism**A. Results**

Examine the six tubes of media, looking for acid and gas. In the presence of acid, bromthymol blue turns yellow. Record your results in the table below. Consult other students for their results and complete the table.

INDIVIDUAL ORGANISMS	LACTOSE		SUCROSE		COMBINATIONS	LACTOSE		SUCROSE	
	Acid	Gas	Acid	Gas		Acid	Gas	Acid	Gas
<i>E. coli</i>					<i>E. coli</i> and <i>P. vulgaris</i>				
<i>P. vulgaris</i>					<i>E. coli</i> and <i>S. aureus</i>				
<i>S. aureus</i>					<i>S. aureus</i> and <i>P. vulgaris</i>				

Bacterial Synergism and Microbial Antagonism

B. Questions

1. Did any of the three organisms produce gas in either lactose or sucrose broth when alone in the medium?

2. Which organisms act synergistically to produce gas in
lactose? _____
sucrose? _____

Ex. 33 Microbial Antagonism

A. Questions

1. Which organisms are antagonistic to *E. coli*? _____

2. Which organisms are antagonistic to *S. aureus*? _____

3. Name an antibiotic substance that is derived from each of the following types of organisms. Also, indicate whether the substance is effective against gram-positive or gram-negative organisms.
A bacterium: _____
An actinomyces: _____
A fungus: _____
4. What role does microbial antagonism play in nature? _____

5. In what physiological way does penicillin affect penicillin-sensitive bacteria? _____

6. How do sulfonamides inhibit bacterial growth? _____

LABORATORY REPORT

34


Student: _____

Desk No.: _____ Section: _____

Temperature: Effects on Growth

A. Pigment Formation and Temperature

1. Draw the appearance of the growth of *Serratia marcescens* on the nutrient agar slants using colored pencils.
2. Which temperature seems to be closest to the optimum temperature for pigment formation? _____
3. What are the cellular substances that control pigment formation and are regulated by temperature? _____



B. Growth Rate and Temperature


If a spectrophotometer is available, dispense the cultures into labeled cuvettes and determine the percent transmittance of each culture. Calculate the O.D. values from the percent transmittance, using the formula given in Exercise 23.

If no spectrophotometer is available, record only the visual reading as +, + +, + + +, and none.

Temp. °C	SERRATIA MARCESCENS			ESCHERICHIA COLI		
	Visual Reading	Spectrophotometer		Visual Reading	Spectrophotometer	
		%T	O.D.		%T	O.D.
5						
25						
38						
42						
55						

Temperature: Effects on Growth

Growth curves of *Serratia marcescens* and *Escherichia coli* as related to temperature.



1. On the basis of the above graph, estimate the optimum growth temperature of the two organisms.

Serratia marcescens: _____

Escherichia coli: _____

2. To get more precise results for the above graph, what would you do?

3. Differentiate between the following:

Thermophile: _____

Mesophile: _____

Psychrophile: _____

4. What is the optimum growth temperature range for most psychrophiles? _____

LABORATORY REPORT

35

Student: _____

Desk No.: _____ Section: _____

Temperature:
Lethal Effects

A. Tabulation of Results

Examine your five Petri plates, looking for evidence of growth. Record on the chalkboard, using a chart similar to the one below, the presence or absence of growth as (+) or (−). When all members of the class have recorded their results, complete this chart.

ORGANISM	60° C					70° C					80° C					90° C					100° C				
	C*	10	20	30	40	C*	10	20	30	40	C*	10	20	30	40	C*	10	20	30	40	C*	10	20	30	40
<i>S. aureus</i>																									
<i>E. coli</i>																									
<i>B. megaterium</i>																									

*C = control tube

1. If they can be determined from the above information, record the **thermal death point** for each of the organisms.

S. aureus: _____ *E. coli*: _____ *B. megaterium*: _____

2. From the following table, determine the **thermal death time** for each organism at the tabulated temperatures.

ORGANISM	THERMAL DEATH TIME				
	60° C	70° C	80° C	90° C	100° C
<i>S. aureus</i>					
<i>E. coli</i>					
<i>B. megaterium</i>					

B. Questions

1. Give three reasons why endospores are much more resistant to heat than are vegetative cells.

- a. _____
- b. _____
- c. _____

Temperature: Lethal Effects

2. Differentiate between the following:

Thermoduric: _____

Thermophilic: _____

3. List four diseases caused by spore-forming bacteria.

a. _____ b. _____
c. _____ d. _____

4. Since boiling water is unreliable in destroying endospores, how should one use heat in medical applications to ensure spore destruction? (three ways)

a. _____
b. _____
c. _____

LABORATORY REPORT

36

Student: _____

Desk No.: _____ Section: _____

pH and Microbial Growth

A. Tabulation of Results

If a spectrophotometer is available, dispense the cultures into labeled cuvettes and determine the percent transmittance of each culture. Calculate the O.D. values from the percent transmittance, using the formula given in Exercise 23. To complete the tables, get the results of the other three organisms from other members of the class, and delete the substitution organisms in the tables that were not used.

If no spectrophotometer is available, record only the visual reading as +, + +, + + +, and none.


pH	<i>Escherichia coli</i>			<i>Staphylococcus aureus</i>		
	Visual Reading	Spectrophotometer		Visual Reading	Spectrophotometer	
		%T	O.D.		%T	O.D.
3						
5						
7						
8						
9						
10						

pH	<i>Alcaligenes faecalis</i> or <i>Sporosarcina ureae</i>			<i>Saccharomyces cerevisiae</i> or <i>Candida glabrata</i>		
	Visual Reading	Spectrophotometer		Visual Reading	Spectrophotometer	
		%T	O.D.		%T	O.D.
3						
5						
7						
8						
9						
10						

pH and Microbial Growth

B. Growth Curves

Once you have computed all the O.D. values on the two tables, plot them on the following graph. Use different colored lines for each species.



C. Questions

1. Which organism seems to grow best in acid media? _____
2. Which organism seems to grow best in alkaline media? _____
3. Which organism seems to tolerate the broadest pH range? _____

LABORATORY REPORT

37, 38

Student: _____

Desk No.: _____ Section: _____

Ex. 37 Osmotic Pressure and Bacterial Growth

A. Results

Record the amount of growth of each organism at the different salt concentrations, using +, ++, +++, and none to indicate degree of growth.

ORGANISM	SODIUM CHLORIDE CONCENTRATION							
	0.5%		5%		10%		15%	
	48 hr	96 hr	48 hr	96 hr	48 hr	96 hr	48 hr	96 hr
<i>Escherichia coli</i>								
<i>Staphylococcus aureus</i>								
<i>Halobacterium salinarium</i>								

B. Questions

1. Evaluate the salt tolerance of the above organisms.

Tolerates very little salt: _____

Tolerates a broad range of salt concentration: _____

Grows only in the presence of high salt concentration: _____

2. How would you classify *Halobacterium salinarium* as to salt needs? Check one.

Obligate halophile _____ Facultative halophile _____

3. Differentiate between the following:

Halophile: _____

Osmophile: _____

4. Supply the following information concerning **mannitol salt agar** (Refer to the Difco Manual).

Composition: _____

For what organism is this medium selective? _____

What ingredient makes it selective? _____

Oligodynamic Action

Ex. 38 Oligodynamic Action

A. Tabulation of Results

Measure the zone of inhibition from the edge of each piece of metal with a millimeter ruler. Record the measurements in the table. Spaces are provided for write-in of additional metals.

METAL	MILLIMETERS OF INHIBITION	
	<i>E. coli</i>	<i>S. aureus</i>
Copper		
Silver		
Aluminum		

B. Questions

1. Which metal seems to exhibit the greatest amount of oligodynamic action?

2. Which metal or metals seem to be ineffective? _____
3. Do these two organisms seem to differ in their susceptibility to oligodynamic action?

Explain: _____

4. What specific chemical substances in bacterial cells are inactivated by heavy metals, affecting growth?

LABORATORY REPORT

39

Student: _____

Desk No.: _____ Section: _____

Ultraviolet Light:
Lethal Effects

A. Tabulation of Results

Your instructor will construct a table similar to the one below on the chalkboard for you to record your results. If substantial growth is present in the exposed area, record your results as + + +. If three or fewer colonies survived, record +. Moderate survival should be indicated as ++. No growth should be recorded as -. Record all information in the table.

ORGANISMS	EXPOSURE TIMES							
<i>S. aureus</i>	10 sec	20 sec	40 sec	80 sec	2.5 min	5 min	10 min	• 20 min
Survival								
<i>B. megaterium</i>	1 min	2 min	4 min	8 min	15 min	30 min	60 min	• 6 min
Survival								

*Plates covered during exposure.

B. Questions

1. What length of time is required for the destruction of non-spore-forming bacteria such as *Staphylococcus aureus*?

2. Can you express, quantitatively, how much more resistant *B. megaterium* spores are to ultraviolet light than *S. aureus* vegetative cells (i.e., how many *times* more resistant are they)?

3. Why is it desirable to remove the cover from the Petri dish when making exposures?

4. In what specific way does ultraviolet light destroy microorganisms?

5. What adverse effect can result from overexposure of human tissues to ultraviolet light?

Ultraviolet Light: Lethal Effects

6. What wavelength of ultraviolet is most germicidal? _____

7. List several practical applications of ultraviolet light to microbial control.

LABORATORY REPORT

40

Student: _____

Desk No.: _____ Section: _____

Evaluation of Disinfectants: The Use-Dilution Method

A. Tabulation of Results

The instructor will draw a table on the chalkboard similar to the one below. Examine your tubes of nutrient broth and pins by shaking them and looking for growth (turbidity). If you are doubtful as to whether growth is present, compare the tubes with a tube of sterile nutrient broth. Record on the chalkboard a plus (+) sign if growth is present and a minus (−) sign if no growth is visible. After all students have recorded their results, complete the following chart.

DISINFECTANT		MINUTES											
		<i>Staphylococcus aureus</i>						<i>Bacillus megaterium</i>					
		C*	1	5	10	30	60	C*	1	5	10	30	60
1:750 Zephiran	Substitution												
5% phenol													
8% formaldehyde													

*C = control tube

B. Questions

1. What conclusions can be drawn from this experiment?

2. Distinguish between the following:

Disinfectant: _____

Antiseptic: _____

Evaluation of Disinfectants: The Use-Dilution Method

3. What factors other than time influence the action of a chemical agent on bacteria?

4. Fill in the equation that explains how the **phenol coefficient** is determined:

P.C. = _____

5. What are some drawbacks that one encounters when attempting to apply the phenol coefficient to all disinfectants?

Student: _____

41

Desk No.: _____ **Section:** _____

Evaluation of Alcohol: Its Effectiveness as a Skin Degerming Agent

A. Tabulation of Results

Count the number of colonies that appear on each of the thumbprints and record them in the following table. If the number of colonies has increased in the second press, record a 0 in percent reduction. Calculate the percentages of reduction and record these data in the appropriate column. Use this formula:

$$\text{Percent Reduction} = \frac{(\text{Colony Count 1st press}) - (\text{Colony Count 2nd press})}{(\text{Colony count 1st press})} \times 100$$

[illegible]

Evaluation of Alcohol: Its Effectiveness as a Skin Degerming Agent

B. Questions

1. In general, what effect does alcohol have on the level of skin contaminants? _____

2. Is there any difference between the effects of dipping versus swabbing? _____
Which method appears to be more effective? _____
3. There is definitely survival of some microorganisms even after alcohol treatment. Without staining or microscopic scrutiny, predict what types of microbes are growing on the medium where you made the right thumb impression after treatment. _____

LABORATORY REPORT

42

Student: _____

Desk No.: _____ Section: _____

Evaluation of Antiseptics: The Filter Paper Disk Method

A. Tabulation of Results

With a millimeter scale, measure the zones of inhibition between the edge of the filter paper disk and the organisms. Record this information. Exchange your plates with other students' plates to complete the measurements for all chemical agents.

DISINFECTANT	MILLIMETERS OF INHIBITION	
	<i>Staphylococcus aureus</i>	<i>Pseudomonas aeruginosa</i>
5% phenol		
5% formaldehyde		
5% iodine		

B. Questions

1. What conclusions can be derived from these results? _____

2. What factors influence the size of the zone of inhibition? _____

Benson: Microbiological Applications Lab Manual, Eighth Edition	Back Matter	Laboratory Reports		© The McGraw-Hill Companies, 2001
--	--------------------	---------------------------	--	--------------------------------------

LABORATORY REPORT

43

Student: _____

Desk No.: _____ Section: _____

Antimicrobial Sensitivity Testing: The Kirby-Bauer Method

A. Tabulation

List the antimicrobics that were used for each organism. Consult tables 43.1 and 43.2 to identify the various disks. After measuring and recording the zone diameters, consult table VII in Appendix A for interpretation. Record the degrees of sensitivity (R, I, or S) in the sensitivity column. Exchange data with other class members to complete the entire chart.

	ANTIMICROBIC	ZONE DIA.	RATING (R, I, S)	ANTIMICROBIC	ZONE DIA.	RATING (R, I, S)
<i>S. aureus</i>						
<i>P. aeruginosa</i>						
<i>Proteus vulgaris</i>						
<i>E. coli</i>						

Antimicrobial Sensitivity Testing: The Kirby-Bauer Method

B. Questions

1. Which antimicrobics would be suitable for the control of the following organisms?
S. aureus: _____
E. coli: _____
P. vulgaris: _____
P. aeruginosa: _____
2. Differentiate between the following:
Narrow spectrum antibiotic: _____

Broad spectrum antibiotic: _____

3. Which antimicrobics used in this experiment would qualify as being excellent broad spectrum antimicrobics?

4. Differentiate between the following:
Antibiotic: _____

Antimicrobial: _____

5. How can drug resistance in microorganisms be circumvented? _____

LABORATORY REPORT

44

Student: _____

Desk No.: _____ Section: _____

Effectiveness of Hand Scrubbing

A. Tabulation of Results

The instructor will draw a table on the chalkboard similar to the one below. Examine the six plates that your group inoculated from the basin of water. Select the two plates of a specific dilution that have approximately 30 to 300 colonies and count all of the colonies of each plate with a Quebec colony counter. Record the counts for each plate and their averages on the chalkboard. Once all the groups have recorded their counts, record the dilution factors for each group in the proper column. To calculate the organisms per milliliter multiply the average count by the dilution factor.


GROUP	0.1 ml COUNT		0.2 ml COUNT		0.4 ml COUNT		DILUTION FACTOR*	ORGANISMS PER MILLILITER
	Per Plate	Average	Per Plate	Average	Per Plate	Average		
A								
B								
C								
D								
E								

*Dilution factors: 0.1 ml = 10; 0.2 ml = 5; 0.4 ml = 2.5

Effectiveness of Hand Scrubbing

B. Graph

After you have completed this tabulation, plot the number of organisms per milliliter that were present in each basin.



C. Questions

1. What conclusions can be derived from this exercise?

2. What might be an explanation of a higher count in Basin D than in B, ruling out contamination or faulty techniques?

3. Why is it so important that surgeons scrub their hands prior to surgery even though they wear rubber gloves?

LABORATORY REPORT

45

Student: _____

Desk No.: _____ Section: _____

Preparation and Care of Stock Cultures

1. Why shouldn't cultures be stored at room temperature or in the incubator for any length of time?

2. Why should stock cultures be reinoculated to new media ("rotated") even if they are stored in the refrigerator? _____

3. For what types of inoculations do you use your
reserve stock culture? _____

working stock culture? _____

4. What is **lyophilization**? _____

What advantage does this procedure have over the method we are using for maintaining stock cultures?

Benson: Microbiological Applications Lab Manual, Eighth Edition	Back Matter	Laboratory Reports		© The McGraw-Hill Companies, 2001
--	--------------------	---------------------------	--	--------------------------------------

LABORATORY REPORT

48–50

Student: _____

Desk No.: _____ Section: _____

Physiological Characteristics of Bacteria

A. Media

List the media that are used for the following tests:

1. Butanediol production
2. Hydrogen sulfide production
3. Indole production
4. Starch hydrolysis
5. Urease production
6. Citrate utilization
7. Fat hydrolysis
8. Casein hydrolysis
9. Catalase production
10. Mixed acid fermentation
11. Glucose fermentation
12. Nitrate reduction

B. Reagents

Select the reagents that are used for the following tests:

- | | |
|-------------------------|---------------------|
| 1. Indole test | Barritt's reagent—1 |
| 2. Voges-Proskauer test | Gram's iodine—2 |
| 3. Catalase test | Hydrogen peroxide—3 |
| 4. Starch hydrolysis | Kovacs' reagent—4 |
| | None of these—5 |

C. Ingredients

Select the ingredients of the reagents for the following tests. Consult Appendix B. More than one ingredient may be present in a particular reagent.

- | | |
|-------------------------|--|
| 1. Oxidase test | α -naphthol—1 |
| 2. Voges-Proskauer test | Dimethyl- α -naphthylamine—2 |
| 3. Indole test | Dimethyl- p -phenylenediamine
hydrochloride—3 |
| 4. Nitrite test | p -dimethylamine benzaldehyde—4 |
| | Potassium hydroxide—5 |
| | Sulfanilic acid—6 |

D. Enzymes

What enzymes are involved in the following reactions?

1. Urea hydrolysis
2. Hydrogen gas production from formic acid
3. Casein hydrolysis
4. Indole production
5. Nitrate reduction
6. Starch hydrolysis
7. Fat hydrolysis
8. Gelatin hydrolysis (Ex. 47)
9. Hydrogen sulfide production

Answers

Media

1. _____
2. _____
3. _____
4. _____
5. _____
6. _____
7. _____
8. _____
9. _____
10. _____
11. _____
12. _____

Reagents

Ingredients

- | | |
|----------|----------|
| 1. _____ | 1. _____ |
| 2. _____ | 2. _____ |
| 3. _____ | 3. _____ |
| 4. _____ | 4. _____ |

Enzymes

1. _____
2. _____
3. _____
4. _____
5. _____
6. _____
7. _____
8. _____
9. _____

Physiological Characteristics of Bacteria

E. Test Results

Indicate the appearance of the following positive test results.

1. Glucose fermentation, no gas
2. Citrate utilization
3. Urease production
4. Indole production
5. Acetoin production
6. Hydrogen sulfide production
7. Coagulation of milk
8. Peptonization in milk
9. Litmus reduction in milk
10. Nitrate reduction
11. Catalase production
12. Casein hydrolysis
13. Fat hydrolysis

Answers

1. _____
2. _____
3. _____
4. _____
5. _____
6. _____
7. _____
8. _____
9. _____
10. _____
11. _____
12. _____
13. _____

F. General Questions

1. Differentiate between the following:

Respiration: _____

Fermentation: _____

Oxidation: _____

Reduction: _____

Catalase: _____

Peroxidase: _____

2. List two or three difficulties one encounters in trying to differentiate bacteria on the basis of physiological characteristics. _____

3. Now that you have determined the morphological, cultural, and physiological characteristics of your unknown, what other kinds of tests might you perform on the organism to assist in identification?

LABORATORY REPORT

52

Student: _____

Desk No.: _____ Section: _____

Enterobacteriaceae Identification: The API 20E System

A. Tabulation of Results

By referring to charts I and II, Appendix D, determine the results of each test and record these results as positive (+) or negative (–) in the table below. Note that the results of the oxidase test must be recorded in the last column on the right side of the table.

ONPG	ADH	LDC	ODC	CIT	H ₂ S	URE	TDA	IND	VP	GEL	GLU	MAN	INO	SOR	RHA	SAC	MEL	AMY	ARA	OXI
1	2	4	1	2	4	1	2	4	1	2	4	1	2	4	1	2	4	1	2	4


NO ₂	N ₂ GAS	MOT	MAC	OF	O	OF	F
1	2	4	1	2	4		

Additional Digits

B. Construction of Seven-Digit Profile

Note in the above table that each test has a value of 1, 2, or 4. To compute the seven-digit profile for your unknown, total up the positive values for each group.

Example:
5 144 572 = *E. coli*



Enterobacteriaceae Identification: The API 20E System

C. Using the API 20E Analytical Index or the API Characterization Chart

If the *API 20E Analytical Index* is available on the demonstration table, use it to identify your unknown, using the seven-digit profile number that has been computed. If no *Analytical Index* is available, use Characterization Chart III in Appendix D.

Name of Unknown: _____

D. Additional Tabulation Blank

If you need another form, use the one below:

API 20E		Reference Number _____		Patient _____		Date _____															
		Source/Site _____		Physician _____		Dept./Service _____															
	ONPG 1	ADH 2	LDC 4	ODC 1	CIT 2	H ₂ S 4	URE 1	TDA 2	IND 4	VP 1	GEL 2	GLU 4	MAN 1	INO 2	SOR 4	RHA 1	SAC 2	MEL 4	AMY 1	ARA 2	OXI 4
5 h																					
24 h																					
48 h																					
Profile Number																					
	NO ₂ 1	GAS 2	MOT 4	MAC 1	OF-O 2	OF-F 4															
5 h																					
24 h																					
48 h																					
Additional Digits																					
Additional Information																					
Identification																					
00-42-012 E-3 (7/80)																					

E. Questions

- What is the intended function of the API 20E system? _____

- In the “real world” who would use this system? _____

- What might be an explanation for the failure of this system to work with some of the bacterial cultures we use? _____

LABORATORY REPORT

53

Student: _____

Desk No.: _____ Section: _____

Enterobacteriaceae Identification: The Enterotube II System

A. Tabulation of Results

Record the results of each test in the following table with a plus (+) or minus (-).

GLUCOSE	GAS PRODUCTION	LYSINE	ORNITHINE	H ₂ S	INDOLE	ADONITOL	LACTOSE	ARABINOSE	SORBITOL	VOGES-PROSKAUER	DULCITOL	PHENYLALANINE DEAMINASE	UREA	CITRATE

B. Identification by Chart Method

If no *Interpretation Guide* is available, apply the above results to chart IV, Appendix D, to find the name of your unknown. Note that the spacing of the above table matches the size of the spaces on chart IV. If this page is removed from the manual, folded, and placed on chart IV, the results on the above table can be moved down the chart to make a quick comparison of your results with the expected results for each organism.

C. Using the Enterotube II Interpretation Guide

If the *Interpretation Guide* is available, determine the five-digit code number by circling the numbers (4, 2, or 1) under each test that is positive, and then totaling these numbers within each group to form a digit for that group. Note that there are two tally charts on the next page of this Laboratory Report for your use.

G L U.	G A S.	L Y S.	O R N.	H ₂ S	I N D.	A D O N.	L A C.	A R A B.	S O R B.	D U L.	P A.	U R E A	C I T.
(2) + (1)	(4) + 2 + 1	4 + (2) + (1)	(4) + (2) + 1	4 + (2) + (1)									
↓	↓	↓	↓	↓									
ID Value	3	4	3	6	3								

The "ID Value" 34363 can be found by thumbing the pages of the *Interpretation Guide*. The listing is as follows:

ID Value
34363

Organism
Klebsiella pneumoniae

Atypical Test Results
None

Conclusion: Organism was correctly identified as *Klebsiella pneumoniae*. In this case, the identification was made independent of the V-P test.

Enterobacteriaceae Identification: The Enterotube II System

D. Tally Charts

ENTEROTUBE® II*

	GLU	GAS	LYS	ORN	H ₂ S	IND	ADON	LAC	ARAB	SORB	DUL	P.A.	UREA	CIT	Confirmatory Test	Result
	2+1		4+2+1			4+2+1		4+2+1		4+2+1		4+2+1				
ID Value	↓ Hexagon ↓ Hexagon ↓ Hexagon ↓ Hexagon ↓ Hexagon															

Culture Number, Case Number or Patient Name

Date

Organism Identified

*VP utilized as confirmatory test only.

ENTEROTUBE® II*

	GLU	GAS	LYS	ORN	H ₂ S	IND	ADON	LAC	ARAB	SORB	DUL	P.A.	UREA	CIT	Confirmatory Test	Result
	2+1		4+2+1			4+2+1		4+2+1		4+2+1		4+2+1				
ID Value	↓ Hexagon ↓ Hexagon ↓ Hexagon ↓ Hexagon ↓ Hexagon															

Culture Number, Case Number or Patient Name

Date

Organism Identified

*VP utilized as confirmatory test only.

E. Questions

- What is the intended function of the Enterotube II System? _____

- In the “real world” who would use this system? _____

- What might be an explanation for the failure of this system to work with some of the bacterial cultures we use? _____

LABORATORY REPORT

54

Student: _____

Desk No.: _____ Section: _____

O/F Gram-Negative Rods Identification: The Oxi/Ferm Tube II System


A. Tabulation of Results and Code Determination

Once you have marked the positive reactions on the side of the tube and circled the numbers that are assigned to each of the positive chambers, as indicated in the example below, add the numbers in each bracketed group to get the five-digit code.

The final step is to look up the code number in the *Oxi/Ferm Tube II Biocode Manual* to determine the genus and species. If confirmatory tests are necessary, the manual will tell you which ones to perform.

In the example below, the code number is 32303. If you look up this number in the *Biocode Manual* you will find on page 25 that the organism is *Pseudomonas aeruginosa*.

Use this procedure to identify your unknown by applying your results to the blank diagrams provided.



B. Results Pads

The diagram shows two identical blank Oxi/Ferm Tube II diagrams for recording results. Each tube has 13 chambers labeled: Ana-Gluc, Arginine, Lysine, Lactose, N₂, Sucrose, Indole, Xylose, Aer-Gluc, Maltose, Mannitol, PA, Urea, Citrate, and Oxidase. Below each chamber, there is a space for a number, and these numbers are grouped by brackets with plus signs. Arrows point from each bracket to a box for the sum of the numbers in that group.

O/F Gram-Negative Rods Identification: The Oxi/Ferm Tube II System

C. Questions

1. What is the intended function of the Oxi/Ferm Tube II System? _____

2. In the “real world” who would use this system? _____

3. What might be an explanation for the failure of this system to work with some of the bacterial cultures we use? _____

LABORATORY REPORT

55

Student: _____

Desk No.: _____ Section: _____

Staphylococcus Identification: The API Staph-Ident System

A. Tabulation of Results


By referring to charts V and VI, Appendix D, determine the results of each test, and record these results as positive (+) or negative (−) in the Profile Determination Table below. Note that two more of these tables have been printed on the next page for tabulation of additional organisms.

	PHS 1	URE 2	GLS 4	MNE 1	MAN 2	TRE 4	SAL 1	GLC 2	ARG 4	NGP 1
RESULTS										
PROFILE NUMBER										

GRAM STAIN		COAGULASE		Additional Information	Identification
MORPHOLOGY		CATALASE			

B. Construction of Four-Digit Profile

Note in the above table that each test has a value of 1, 2, or 4. To compute the four-digit profile for your unknown, total up the positive values for each group.



C. Final Determination

Refer to the Staph-Ident Profile Register (chart VII, Appendix D) to find the organism that matches your profile number. Write the name of your unknown in the space below and list any additional tests that are needed for final confirmation. If the materials are available for these tests, perform them.

Name of Unknown: _____

Additional Tests: _____

Styphylococcus Identification: The API Staph-Ident System

	PHS 1	URE 2	GLS 4	MNE 1	MAN 2	TRE 4	SAL 1	GLC 2	ARG 4	NGP 1
RESULTS										
PROFILE NUMBER										
GRAM STAIN				Additional Information				Identification		
MORPHOLOGY										
COAGULASE										
CATALASE										

	PHS 1	URE 2	GLS 4	MNE 1	MAN 2	TRE 4	SAL 1	GLC 2	ARG 4	NGP 1
RESULTS										
PROFILE NUMBER										
GRAM STAIN				Additional Information				Identification		
MORPHOLOGY										
COAGULASE										
CATALASE										

D. Questions

- What is the intended function of the API Staph-Ident System? _____

- In the “real world” who would use this system? _____

- What might be an explanation for the failure of this system to work with some of the bacterial cultures we use? _____

LABORATORY REPORT

56, 57

Student: _____

Desk No.: _____ Section: _____

Ex. 56 Microbial Population Counts of Soil

A. Tabulation of Results

Select the best plate and count the organisms on a Quebec colony counter. Tabulate your results and the results of other students near you who cultured the other types of soil organisms. Be sure to count only representative types.

ORGANISMS	COUNT PER PLATE	DILUTION	ORGANISMS PER GRAM OF SOIL
Bacteria			
Actinomycetes			
Molds			

B. Conclusions

What generalizations can you make from this exercise?

C. Questions

1. Why would the number of bacteria present in the soil actually be higher than the number determined by your plate count? _____

2. What other types of microbes might be present in your soil sample? _____

Isolation of an Antibiotic Producer from Soil

Ex. 57 Isolation of an Antibiotic Producer from Soil

A. Results

Describe in detail how the experiment turned out.

B. Questions

1. What four genera of microbes produce the vast majority of antibiotics? _____

2. What is known about the importance of these antibiotics in the natural habitat of these microbes?

LABORATORY REPORT

59

Student: _____

Desk No.: _____ Section: _____

Nitrogen Fixing Bacteria


A. Questions

1. What enzyme is responsible for nitrogen fixation? _____
2. What metal is essential for this enzyme to function? _____
3. Why is nitrogen fixation so important? _____

B. Azotobacteraceae

Identify the characteristics of your isolate that led to your decision as to its identification.

1. Were the colonies pigmented? _____
If pigmented, what color? _____
Was fluorescence present? _____
If fluorescent, what color? _____
2. Were you able to see cysts? _____
3. Is the organism motile? _____
4. Gram reaction? _____
5. Any other pertinent characteristic? _____



Azotobacteraceae

Name of Organism: _____


Draw a few cells of organisms in circle at right.

C. Rhizobiaceae

Draw a few representative cells of *Rhizobium* in the circle to the right.

1. What is the most important criterion for species identification in the family of nitrogen-fixers?

2. From the standpoint of amount of nitrogen fixation, is this group of nitrogen-fixers **more** or **less important** than the Azotobacteraceae?



Rhizobiaceae

Benson: Microbiological Applications Lab Manual, Eighth Edition	Back Matter	Laboratory Reports		© The McGraw-Hill Companies, 2001
--	--------------------	---------------------------	--	--------------------------------------

LABORATORY REPORT

60

Student: _____

Desk No.: _____ Section: _____

Ammonification in Soil

A. Tabulation of Results

Record the presence or absence of ammonia in the tubes of media:

- + slight ammonia (faint yellow)
- ++ moderate ammonia (deep yellow)
- +++ much ammonia (brown precipitate)
- no ammonification

INCUBATION TIME	AMOUNT OF AMMONIA		pH	
	Control	Peptone with Soil	Control	Peptone with Soil
4 days				
7 days				

B. Questions

1. In the natural situation, what compounds serve as the source of the ammonia released in ammonification?

2. In simple terms, what occurs during decay? _____

3. Differentiate between the following:

Peptonization: _____

Ammonification: _____

4. What happens to ammonia in soil when **denitrification** of ammonia takes place? _____

Benson: Microbiological Applications Lab Manual, Eighth Edition	Back Matter	Laboratory Reports		© The McGraw-Hill Companies, 2001
--	--------------------	---------------------------	--	--------------------------------------

LABORATORY REPORT

61, 62

Student: _____

Desk No.: _____ Section: _____

Ex. 61 Isolation of a Denitrifier from Soil Using Nitrate Agar

A. Observations:

1. **Second Period.** Describe the types of colonies observed growing on the nitrate agar plate. _____

2. **Third Period.** What evidence do you have of the presence of a denitrifier in the soil sample you are working with? _____

3. **Microscopic Examination.** Describe the appearance of organisms that were observed on a gram-stained slide. _____

4. **Further Testing.** If you performed any further tests for species identification, state any conclusions that were made. _____

B. Questions:

1. Why are denitrifying bacteria costly to farmers? _____

2. Why are denitrifying bacteria essential to the existence of life on plant earth? _____

3. Of what value is denitrification to the organism? _____

Isolation of a Denitrifier from Soil: Using an Enrichment Medium

Ex. 62 Isolation of a Denitrifier from Soil: Using an Enrichment Medium

A. Observations:

1. **First Enrichment Culture Appearance.** Describe the appearance of the first enrichment culture in the second laboratory period: _____

2. **Microscopic Examination:** Describe the appearance of the cells as seen on:
Wet Mount: _____

Gram-stained Slide: _____

3. **Second Enrichment Culture Appearance.** Describe the appearance of the second enrichment culture after five days incubation: _____

4. **Agar Plate Colonies.** Describe the characteristics of the colonies on the nitrate succinate-mineral salts agar: _____

Do these characteristics match those of *Paracoccus denitrificans*? _____
5. **Final Gram-stained Slide.** How do the microscopic characteristics of your organism match the description of *Paracoccus denitrificans*? _____

B. Questions: Answer the questions that are provided on the front page of this Laboratory Report.

LABORATORY REPORT

63

Student: _____

Desk No.: _____ Section: _____

Bacteriological Examination of Water: Qualitative Tests

A. Results of Presumptive Test (MPN Determination)

Record the number of positive tubes on the chalkboard and on the following table. When all students have recorded their results with the various water samples, complete this tabulation. Determine the MPN according to the instructions on page 224.

WATER SAMPLE (SOURCE)	NUMBER OF POSITIVE TUBES				MPN
	3 Tubes DSLB 10 ml	3 Tubes SSLB 1.0 ml	3 Tubes SSLB 0.1 ml	3 Tubes SSLB 0.01 ml	

B. Results of Confirmed Test

Record the results of the confirmed tests for each water sample that was positive on the presumptive test.

WATER SAMPLE (SOURCE)	POSITIVE	NEGATIVE